

Recommend Clinical Supplies Size

Item	Adult	PediaSIM	BabySIM	iStan
Urinary Catheter	14-16 Fr	10 Fr	8 Fr	16 Fr
Nasogastric Tube	14 Fr	10 Fr	8 Fr	14 Fr
ETT	7.5-8.0	5	3.5	7.5-8.0
LMA	3	2	1.5	3
Oropharyngeal Airway	90 mm	60 mm	43 mm	90 mm
Nasopharyngeal Airway	30 mm	24 mm	20 mm	30 mm
IV Cannula	20-22 gauge	20-22 gauge	20-22 gauge	20-22 gauge
Chest Tube	28 Fr	20 Fr	12 Fr	28 Fr
Needle Decompression	14 gauge 3-6 cm long	Same	Same	14 gauge 3-6 cm long
Pericardiocentesis	16-18 gauge 10-15 cm long	Same	Same	N/A
Combitube	37 Fr	N/A	N/A	37 Fr